

## Neonatal Weight Loss Nomogram Tool (NEWT) Reference

To use the Neonatal Weight Loss Nomogram Tool (NEWT), the patient must be **less than 30 days in age** AND the following fields **must be populated** in the patient's chart:

- Birth Date & Time
- Delivery Type (vaginal or C-section)
- Feeding Plan for Newborn (breastmilk, formula, or bottle)
- Birth Weight
- At least one subsequent weight

### Accessing the Neonatal Weight Loss Nomogram Tool (NEWT)

Intake and Output  
Labs  
Lines/Tubes/Drains  
Media Gallery  
Microbiology  
MoodTrek Overview  
Neonate Overview  
Neonate Transfusion  
Neonate Weights and Measurements  
New Order Entry  
Newborn Weight Tool (NEWT)  
Order Profile

Components  
Clear Preferences


Add the **NEWT** section by marking it for display within the *Customize Menu (pancake) / Components* on the Workflow page

### Main Features of the Tool:



Changes made in this component will not be saved.

Newborns charted with feeding method of breastmilk and formula will be displayed as being breastfed only.



**newt**

Birth Details

Weight: 3.404 kg    Date: Sep 23    Time: 06:29

Vaginal    [Edit Details](#)

Measurements

Day	Weight	Change	Add New
Birth	3.404 kg	—	<a href="#">Edit</a>
1	3.263 kg	-4.1%	<a href="#">Edit</a> <a href="#">x</a>
1.9	3.138 kg	-7.8%	<a href="#">Edit</a> <a href="#">x</a>
4.2	3.22 kg	-5.4%	<a href="#">Edit</a> <a href="#">x</a>

**PennState Health**  
Children's Hospital

If the patient is between 72 hours and 30 days of age, a 30-day view will display